

# GX Wolf Microscopes


product datasheet

## Quality

Microscope Slides · Coverslips · Prepared Slides

### GXM Glass Slides Specifications

#### Microscope Slides



Item		Specification and Packing	
<b>7101</b> Plain Slide Ground Edges		Size mm	25.4 X 76.2
		Thickness(mm)	1.0-1.2
			0.8-1.0
		Packing (PCS/Box)	50 72
Boxes per Carton	50		
<b>7102</b> Plain Slide Unground Edges		Size mm	25.4 X 76.2
		Thickness(mm)	1.0-1.2
			0.8-1.0
		Packing (PCS/Box)	50 72
Boxes per Carton	50		
<b>7103</b> Single Concave Slide Ground Edges		Size mm	25.4 X 76.2
		Thickness(mm)	1.1-1.3
		Packing (PCS/Box)	50 72
			Boxes per Carton

# GXM Glass Slides Specifications

## Microscope Slides



7101 & 7102



7103



7104



7105



7105-1







Box of 72

Item		Specification and Packing			
<b>7104</b> Double Concave Slide Ground Edges		Size mm	25.4 X 76.2		
		Thickness(mm)	1.1-1.3		
		Packing (PCS/Box)	<table border="1" style="margin: auto;"> <tr><td style="text-align: center;">50</td></tr> <tr><td style="text-align: center;">72</td></tr> </table>	50	72
		50			
72					
Boxes per Carton	50				
<b>7105</b> Markable Slide Frosted One End One Side Ground Edges		Size mm	25.4 X 76.2		
		Thickness(mm)	<table border="1" style="margin: auto;"> <tr><td style="text-align: center;">1.0-1.2</td></tr> <tr><td style="text-align: center;">0.8-1.0</td></tr> </table>	1.0-1.2	0.8-1.0
		1.0-1.2			
		0.8-1.0			
Packing (PCS/Box)	<table border="1" style="margin: auto;"> <tr><td style="text-align: center;">50</td></tr> <tr><td style="text-align: center;">72</td></tr> </table>	50	72		
50					
72					
Boxes per Carton	50				
<b>7105-1</b> Markable Slide Frosted One End On Side Unground Edges		Size mm	25.4 X 76.2		
		Thickness(mm)	<table border="1" style="margin: auto;"> <tr><td style="text-align: center;">1.0-1.2</td></tr> <tr><td style="text-align: center;">0.8-1.0</td></tr> </table>	1.0-1.2	0.8-1.0
		1.0-1.2			
		0.8-1.0			
Packing (PCS/Box)	<table border="1" style="margin: auto;"> <tr><td style="text-align: center;">50</td></tr> <tr><td style="text-align: center;">72</td></tr> </table>	50	72		
50					
72					
Boxes per Carton	50				

# GXM Glass Slides Specifications

## Microscope Slides



Item		Specification and Packing						
<p><b>7106</b> Markable Slide</p> <p>Frosted Both Ends</p> <p>One Side Ground Edges</p>		Size mm	25.4 X 76.2	<table border="1"> <tr> <td data-bbox="857 527 1105 583">Thickness(mm)</td> <td data-bbox="1105 527 1424 583">1.0-1.2</td> </tr> <tr> <td data-bbox="857 583 1105 621"></td> <td data-bbox="1105 583 1424 621">0.8-1.0</td> </tr> </table>	Thickness(mm)	1.0-1.2		0.8-1.0
Thickness(mm)	1.0-1.2							
	0.8-1.0							
<p><b>7107</b> Markable Slide</p> <p>Frosted One End Both Sides</p> <p>Ground Edges</p>		Size mm	25.4 X 76.2	<table border="1"> <tr> <td data-bbox="857 846 1105 903">Thickness(mm)</td> <td data-bbox="1105 846 1424 903">1.0-1.2</td> </tr> <tr> <td data-bbox="857 903 1105 940"></td> <td data-bbox="1105 903 1424 940">0.8-1.0</td> </tr> </table>	Thickness(mm)	1.0-1.2		0.8-1.0
Thickness(mm)	1.0-1.2							
	0.8-1.0							
<p><b>7107-1</b> Markable Slide</p> <p>Frosted One End Both Sides</p> <p>Unground Edges</p>		Size mm	25.4 X 76.2	<table border="1"> <tr> <td data-bbox="857 1171 1105 1228">Thickness(mm)</td> <td data-bbox="1105 1171 1424 1228">1.0-1.2</td> </tr> <tr> <td data-bbox="857 1228 1105 1266"></td> <td data-bbox="1105 1228 1424 1266">0.8-1.0</td> </tr> </table>	Thickness(mm)	1.0-1.2		0.8-1.0
Thickness(mm)	1.0-1.2							
	0.8-1.0							
		Packing (PCS/Box)	50 72	<table border="1"> <tr> <td data-bbox="857 1323 1105 1379">Boxes per Carton</td> <td data-bbox="1105 1323 1424 1379">50</td> </tr> </table>	Boxes per Carton	50		
Boxes per Carton	50							

# GXM Glass Slides Specifications

## Microscope Slides



7101 & 7102



7103




7104



7105



7105-S



Box of 72

Item		Specification and Packing	
<b>7108</b> High-Quality Slide  Frosted Both Ends Both Sides  Ground Edges		Size mm	25.4 X 76.2
		Thickness(mm)	1.0-1.2 0.8-1.0
		Packing (PCS/Box)	50 72
		Boxes per Carton	50
		<b>7109</b> High Quality Slide  Color Frosted  One End One Side  Ground Edges	
Thickness(mm)	1.0-1.2 0.8-1.0		
Packing (PCS/Box)	50 72		
Boxes per Carton	50		
<b>7110</b> Markable Slide  Frosted One side  Ground Edges			
		Thickness(mm)	1.0-1.2 0.8-1.0
		Packing (PCS/Box)	50 72
		Boxes per Carton	50
		<b>7111</b> Markable Slide  Frosted Both Sides  Ground Edges	
Thickness(mm)	1.0-1.2		
Packing (PCS/Box)	50 72		
Boxes per Carton	50		

# GX Wolf Microscopes

product datasheet

## GXM Coverslips Specifications

### Coverslips



### Coverslips

Type	Dimensions mm	Thickness mm	Product Code	Qty per box
Square	18x18	0.13-0.17	18	100
Square	20x20	0.13-0.17	20	100
Square	22x22	0.13-0.17	22	100
Square	24x24	0.13-0.17	24	100
Rectangular	24x32	0.13-0.17	24/32	100
Rectangular	24x40	0.13-0.17	24/40	100
Rectangular	24x50	0.13-0.17	24/50	100
Rectangular	24x60	0.13-0.17	24/60	100
Circular	13	0.13-0.17	13C	100
Circular	15	0.13-0.17	15C	100
Circular	16	0.13-0.17	16C	100
Circular	18	0.13-0.17	18C	100
Circular	19	0.13-0.17	19C	100
Circular	20	0.13-0.17	20C	100

### Wolf Laboratories Limited

Colenso House  
1 Deans Lane  
Pocklington  
York YO42 2PX

Tel 01759 301142  
Fax 01759 301143  
Email [sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)  
Web [www.wolflabs.co.uk](http://www.wolflabs.co.uk)