


GX Wolf Microscopes

product datasheet

Stage Graticules



Transmitted Light Stage Graticule Glass Slides		
Product Code	Scale Size	Scale Interval (mm) (1000 microns per mm)
GXMG 1	10mm	0.1mm
GXMG 2	5mm	0.05mm
GXMG 4	0.1 inch	0.001 inch
GXMG 8	1mm	0.01mm
GXMG 11	0.005 inch	0.0001 inch
GXMG 12	0.1mm	0.002mm
GXMG 16	1mm	0.01mm cross
GXMG 22	2mm	0.01mm

Reflected Light Stage Graticule Slides		
Product Code	Scale Size	Scale Interval (mm) (1000 microns per mm)
GXMGR 1	10mm	0.1mm
GXMGR 2	5mm	0.05mm
GXMGR 4	0.1 inch	0.001 inch
GXMGR 8	1mm	0.01mm
GXMGR 11	0.005 inch	0.0001 inch
GXMGR 12	0.1mm	0.002mm
GXMGR 16	1mm	0.01mm
GXMGR 22	2mm	0.01mm

Wolf Laboratories Limited
Colenso House
1 Deans Lane
Pocklington
York YO42 2PX

Tel 01759 301142
Fax 01759 301143
Email sales@wolflabs.co.uk
Web www.wolflabs.co.uk

Due to our policy of continuous development specifications may change without notice

June 05