

# GX Wolf Microscopes

## product datasheet



### **GXM Ergonomic, Motorised, Height Adjustable, Stable Microscope Workstation Tables**

Wolf Laboratories microscope workstation tables are specifically designed to provide a stable and versatile platform for an optical microscope.

They are available in 2 styles, the GXMAT16 model is equipped with lockable wheels and the GXMCOS560 is a larger corner format. Both have a component drawer and conveniently located controls on the underside of the table top for raising or lowering the work surface height effortlessly with quiet running motors.

These tables are both a more versatile and often a more economic alternative to an ergonomic head on a microscope. In fact the GXMAT16 costs less than most standard workstation tables which have no motorisation.

They are constructed from a tough, resistant melamine-type table top and a metal stand. They can be easily serviced by an electrical engineer.

Model	GXMAT16	GXMCOS560
Elevation Range (mm)	640-800	640-800
Maximum Load	45	45
Worksurface Size (mm)	950x950	1450max X 950

**Wolf Laboratories Limited**  
Colenso House  
1 Deans Lane  
Pocklington  
York YO42 2PX

Tel 01759 301142  
Fax 01759 301143  
Email [sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)  
Web [www.wolflabs.co.uk](http://www.wolflabs.co.uk)

Due to our policy of continuous development specifications may change without notice

April 08