

# GX Wolf Microscopes

## Product datasheet

### Track Stands

Stand and illumination systems for SZ, PZ and MZ series microscopes

GX Microscopes Tracks Stands for Z series Microscopes					
Max Weight	4kg	4kg	4kg	4kg	4kg
Max Size of Base	320x305x26mm	320x305x26mm	320x305x26mm	320x305x26mm	320x305x26mm
Maximum Distance from vertical post to the front of base	273mm	273mm	273mm	273mm	273mm
Diameter of Post	32mm	32mm	32mm	32mm	32mm
Height of Post	325mm	325mm	325mm	325mm	325mm
Type and Diameter of Microscope holder ring/focus mount	For SZ series 76mm	For SZ series 39mm	For SZ series 76mm	For SZ series 76mm Coarse and Fine Focus	For PZ series
Focus Range of Focus Mechanism	25mm	25mm	25mm	25mm	25mm
Illumination	None – requires accessory illumination	None – requires accessory illumination	Top and base illumination	None – requires accessory illumination	None – Requires Accessory Illumination
Power	None	None	85-265V 50/60Hz	None	None

#### Wolf Laboratories Limited

Colenso House  
1 Deans Lane  
Pocklington  
York YO42 2PX

Tel 01759 301142  
Fax 01759 301143  
Email [sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)  
Web [www.wolflabs.co.uk](http://www.wolflabs.co.uk)

Due to our policy of continuous development specifications may change without notice

June 04