

GX Wolf Microscopes


product datasheet

XTD-20 Series Stereo Microscopes


High Quality
Optics, Robust
Microscopes ideal
for Educational
and Hobby use



XTD-20



XTD-20L



XTD-20F

XTD-20 Series Specifications

	XTD-20	XTD-20L	XTD-20F
Eyepieces	WF10X	WF10X	WF10X
Objective	2X WD 81mm	2X WD 81mm	2X WD 81mm
Total Magnification	20X	20X	20X
Stage Plate	Reversible white/black	Reversible white/black	Reversible white/black Frosted Glass
Light Source	None	Top Lamp	Top and Bottom Lamp
Options			
Eyepieces	WF5X/18mm WF5X/13mm WF20X/10mm		
Objectives	1X WD 56mm 3X WD 61mm 4X WD 58mm 6X WD 66mm		
Other	Bulbs, Stage plates, Wooden case (lockable)		

The XTD-20 series of microscopes are superb quality, robust instruments ideal for use in education and for hobbies. They are equipped with first class quality optics normally provided in instruments costing many times more. The result is a crystal clear, bright and perfectly coloured image. The long working distance makes it easy to manipulate specimens under the microscope. The focus controls are small and easy for children to manipulate. The all metal construction results in a long life for the microscope. A range of accessory optics and parts are available to expand the capabilities of each model.

Wolf Laboratories Limited
Colenso House
1 Deans Lane
Pocklington
York YO42 2PX

Tel 01759 301142
Fax 01759 301143
Email sales@wolflabs.co.uk
Web www.wolflabs.co.uk

Due to our policy of continuous development specifications may change without notice

July 05