

Single Channel Manual Pipettors



Discovery+, HTL

- Accurate, precise and robust
- Light and ergonomically designed. NEW designed spring system has reduced plunger pressures by up to 50%
- Fully autoclavable without the need to recalibrate*
- Ultra low pipetting and blowout forces
- Volume can be adjusted either with pushbutton or thumb wheel
- The top of the pushbutton rotates independently preventing accidental volume changing
- The tip ejector is adjustable to a range of tips including filter and gel loading tips
- Compatible with Gilson and a wide range of pipette tips
- Highly durable, colour coded shafts
- Colour coded rings for technique or user identification
- UV resistant
- Supplied with calibration tool and pipettor hanger
- 3 year warranty

*It is recommended that pipettor accuracy is checked after 10 autoclave cycles

HDM020	0.2-2.0ul HTL Discovery+ pipettor
HDM021	0.5-10ul HTL Discovery+ pipettor
HDM022	2.0-20ul HTL Discovery+ pipettor
HDM023	5.0-50ul HTL Discovery+ pipettor
HDM024	10-100ul HTL Discovery+ pipettor
HDM025	20-200ul HTL Discovery+ pipettor
HDM026	50-250ul HTL Discovery+ pipettor
HDM027	100-1000ul HTL Discovery+ pipettor
HDM028	1000-5000ul HTL Discovery+ pipettor
HDM029	1000-10000ul HTL Discovery+ pipettor