

eCount™ Colony Counter

Ideal for rapid, accurate counts

All-in-one eCount combines a multi-function electronic counter with a Sharpie® pen for marking Petri dishes to prevent missing or double-counting colonies. Light, ergonomic body is balanced for easy manipulation.



Touch controls for multiple counts, sums, counting up or down, and more! Counter automatically beeps with each count—or select the LED flash mode for a visual alert of each count. Press the lock key to write with Sharpie pen without counting—makes writing count totals simple.

Sharpie pen provides highly visible marks on both plastic and glass dishes. Counter accepts both fine and extra fine point Sharpie pens. Easily change colors for separate counts. Includes one fine-point Sharpie pen and weighted stand.



Cat No.	Description	Units
HS15092A	Electronic Colony Counter w/Sharpie Pen	Each
HS15092ZP	Replacement Pen for eCount	Pckg. of 12