

WOLF LABORATORIES

Tel: 01759 301142 Fax: 01759 301143 Email: sales@wolflabs.co.uk
Website: www.wolflabs.co.uk

Foam Culture Tube Rack



Floats in water bath

This culture tube rack is constructed of chemical and solvent resistant polyester foam. Tubes fit securely in the sockets.
Packs of 6.

Cat No.	Wells	Tubes	Units
ZXZHS21645A	50	10x75, 12x75 or 13x100mm	Pck. of 6
ZXZHS21645C	30	250 and 400ul	Pck of 6