

# WOLF LABORATORIES

Tel: 01759 301142 Fax: 01759 301143 Email: [sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)

Website: [www.wolflabs.co.uk](http://www.wolflabs.co.uk)

## Economy Microscope Slide Boxes



Constructed of durable polypropylene, slide boxes hold 25 or 100 standard slides (75 x 25 mm and 3 x 1").

Hinged lids are held securely closed by latch lock. Stacking feet maximize stability No liner.

The 25-place box measures 141 x 92 x 36 mm.

The 100-place box measures 210 x 169 x 37 mm.

<b>Cat No</b>	<b>Capacity</b>	<b>Colour</b>
<b>HS15990A</b>	25 Place	Blue
<b>HS15990B</b>	25 Place	Red
<b>HS15990C</b>	25 Place	White
<b>HS15991C</b>	100 Place	White