

WOLF LABORATORIES

Tel: 01759 301142 Fax: 01759 301143 Email: sales@wolflabs.co.uk
Website: www.wolflabs.co.uk

Microscope Slide Dispenser



Unique design makes it easy to load and dispense slides.
Standard slides (75 x 25 mm and 3 x 1") fit effortlessly.
Turn either side wheel to dispense slides one at a time.
Polystyrene dispenser holds up to 91 slides.
Measures 113 x 123 x 163 mm.
Product Code ZXZHS15984