

# WOLF LABORATORIES

Tel: 01759 301142 Fax: 01759 301143 Email: sales@wolflabs.co.uk  
Website: www.wolflabs.co.uk

## Microscope Slide Mailers

### 2 Place Slide Mailers



These 2-place slide protectors hold standard (75 x 25 mm and 3 x 1") slides securely. With large, easy-to-secure snaps, our polypropylene mailers provide safe storage for valuable slides. Array of colours brings some fun to the lab and is useful for colour coding. Assorted colours: red, yellow, blue, purple, and orange. Measures 93 x 68 x 5 mm. Pack of 25.

Cat No.	Description	Colour	Units
<b>ZXZHS15983A</b>	2-Place Slide Mailer	Assorted	Pk of 25
<b>HS15983B</b>	2-Place Slide Mailer	Natural	Pk of 25

### 5 Place Slide Mailers



Ideal for shipping, storing, or staining slides. These polypropylene containers hold five standard slides (75 x 25 mm and 3 x 1"). Inside wall features slotted channels to keep slides separated and to make removal easier. Close securely to keep slides safe. The side-opening mailer measures 81 x 19 x 31 mm; end-opening mailer is 43 x 24 x 88 mm. Pack of 25.

Cat No.	Description	Colour
<b>HS15982</b>	Side-Opening Mailer	Natural
<b>HS15986</b>	End-Opening Mailer	Natural