

WOLF LABORATORIES

Tel: 01759 301142 Fax: 01759 301143 Email: sales@wolflabs.co.uk
Website: www.wolflabs.co.uk

Disposable Slide Mailers



1 or 2-place mailers

Heavy cardboard holders with thumb grooves make the removal of slides easy. We offer 1-slide and 2-slide versions. Each will hold standard (3 x 1", or 75 x 25mm) slides. The one-place mailer measures 95 x 40 x 5mm. The two-place version is 103 x 79 x 5mm.

Cat No.	Description
---------	-------------

ZXZHS9903	1-Place Cardboard Slide Mailer pk/25
-----------	--------------------------------------