

WOLF LABORATORIES

Tel: 01759 301142 Fax: 01759 301143 Email: sales@wolflabs.co.uk
Website: www.wolflabs.co.uk

Economy Counters



Mechanical, benchtop counters are ideal for determining the number of cells, counting parts, etc. Use with blood cell label or insert your own. Each key records up to 999. Total window keeps track of the total number of stokes made on all other keys; at one hundred totals, counter rings bell. Counter with 6-key measures 250 x 76 x 56mm. 9-key counter measures 320 x 76 x 56mm.

Cat No.	Description	Units
ZXZHS5360A	6-Key Economy Counter	Each
ZXZHS5360B	9-Key Economy Counter	Each