


# WOLF LABORATORIES

Tel: 01759 301142 Fax: 01759 301143 Email: sales@wolflabs.co.uk  
Website: www.wolflabs.co.uk

## PCR® Transfer Tray



### Handy when using with thermal cyclers

A time-saver for anyone working with individual 0.2 mL tubes or strips of 8 or 12 tubes. Moving the 0.2 mL tubes to the thermal cycler is super-simple with this tray—just use handles to lift and carry. Transfer tray is a thin strip of plastic designed to rest perfectly over a 96-well tube rack. Made out of high-impact polypropylene. Available in natural. Tray measures 114 x 74 x 11 mm. Pack of 5.

Cat No.	Description	Description
ZXZHS23465B	Natural	Pack of 5