

WOLF LABORATORIES

Tel: 01759 301142 Fax: 01759 301143 Email: sales@wolflabs.co.uk
Website: www.wolflabs.co.uk

Pocket Timer/Stopwatch



Four independent timing channels with visual and audible alarms. Simply move the program switch to use as a stopwatch, a count down timer (up to 100 hours), or a clock. Read out time on 6-digit LCD, 0.25"H. Use spring clip to attach to lab coats, magnetic back to put on metal cabinets, or flip open the stand to place on benchtop. Includes two-year silver oxide battery.
Measures 71 x 61 x 25mm

Cat No.	Description	Units
ZXZHS24490D	Pocket Timer/Stopwatch	Each