

# WOLF LABORATORIES

Tel: 01759 301142 Fax: 01759 301143 Email: [sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)  
Website: [www.wolflabs.co.uk](http://www.wolflabs.co.uk)

## Rack for Conical Tube – Snap Together



Easy-to-assemble polypropylene racks snap together securely  
Designed for stable use in water baths  
Features alphanumeric grid reference aids for tube identification  
Autoclavable  
Shipped flat  
Wells for 15 mL tubes measure 17 mm dia  
Wells for 50 mL tubes measure 28 mm dia

Cat No.	Tubes	Wells	Colour	L x W x D"
ZXZHS24306	15mL & 50mL	30 & 20	Orange	207 x 171 x 60