

# WOLF LABORATORIES

Tel: 01759 301142 Fax: 01759 301143 Email: sales@wolflabs.co.uk  
Website: www.wolflabs.co.uk

## Reversible PCR® Rack



**Removable hinged lid fits both sides**

Flip this polypropylene reversible rack to the size of tubes you need. PCR side of the rack has 168 wells that hold 0.2 mL tubes. The wells accommodate 8- or 12-tube strips or individual tubes. The opposite side of the rack has 40 wells that hold 0.5 mL tubes. Both sides of the rack have 12 wells for 1.5 mL tubes. Autoclavable. Available in assorted colours: blue, green, purple, yellow, and orange.

Rack measures 206 x 131 x 54 mm.

Part Number	Description	Units
ZXZHS2344A	Reversible PCR Rack, Assorted Colours	Pck of 5