

# WOLF LABORATORIES

Tel: 01759 301142 Fax: 01759 301143 Email: sales@wolflabs.co.uk  
Website: www.wolflabs.co.uk

## Reversible Tube Rack 96 Well



### Easy-to-grasp tabs

Polypropylene rack has 96 wells for 0.5 ml microcentrifuge tubes on one side and 96 wells for 1.5 to 2.0 ml microcentrifuge tubes on the other side

Wells are arranged in standard 8 x 12 array

Handy lid snaps into place via tabs on each end of rack

0.5 ml wells are separated by moulded-in grid lines, making them easier to see, and raised rings are around the 1.5 to 2.0 ml wells make tubes easier to remove

Autoclavable

Rack measures 246 x 121 x 50 mm

Cat No.	Description	Units
HS2345B	96-Well Reversible Racks, Natural	Pck. of 5