

WOLF LABORATORIES

Tel: 01759 301142 Fax: 01759 301143 Email: sales@wolflabs.co.uk
Website: www.wolflabs.co.uk

Solution Basins



Designed and packaged for superior handling!

Our redesigned reservoirs are improved in every way _ with thicker walls, crush-resistant "turrets" at the corners, improved pour spouts, 25 mL graduation marks, and tabs to prevent stacked reservoirs from sticking together. The bottom of each reservoir tapers down into a narrow trough, to help you get the last drop of liquid.

Sterile reservoirs are sold in packs of 50. Individually sealed in our new, Sani-Tear, high density polyethylene bag and gamma-irradiated in an ISO-certified US facility.

Non-sterile reservoirs are sold in packs of 100

Cat No.	Volume	Material	Sterility
HS20521A	55 mL	Clear PVC	Non-Sterile
HS20521B	55 mL	White Polystyrene	Non-Sterile
HS20621C	55 mL	White Polystyrene	Sterile
HS20721B	100 mL	White Polystyrene	Non-Sterile
HS20721C	100 mL	White Polystyrene	Sterile