

WOLF LABORATORIES

Tel: 01759 301142 Fax: 01759 301143 Email: sales@wolflabs.co.uk
Website: www.wolflabs.co.uk

Test Tube Racks - 1 Well



ZXZHS23051



ZXZHS23053A



ZXZHS23052

Keep single tubes upright in confined spaces

Select among our three one-well racks--epoxy-coated steel wire racks, friction-fit polystyrene racks, or interlocking polypropylene racks.

ZXZHS23051

The autoclavable, epoxy-coated steel wire rack provides more friction with the benchtop for more secure positioning. Both racks hold conical- or round-bottom tubes. Friction-fit rack measures 92 mm dia x 30 mm H. Coated wire rack measures 90 x 67 x 67 mm. Friction-fit racks sold in pack of 5.

ZXZHS23053A

The friction-fit polystyrene rack slides smoothly across the benchtop and snugly fits tubes from 27.5 to 29.0 mm in diameter.

ZXZHS23052

Interlock polypropylene racks to form a line of 50 mL tubes, or use as a single, stand-alone rack. Molded circular base ring keeps conical or round-bottom tubes in place. Autoclavable. Assorted colors: one each of blue, green, yellow, pink, and orange. Measure 83 x 70 x 67 mm. Pack of 5.

Cat No.	Description	Tubes	Wells	Color	Units
ZXZHS23051	Epoxy-Coated Wire Rack	50 mL	1	Yellow	Each
ZXZHS23053A	Polystyrene Friction-Fit Rack	50 mL	1	Blue	Pk. of 5
ZXZHS23052	1 - Well 50 mL Tube Racks	50 mL	1 x 3	Assorted	Pk. of 5