

WOLF LABORATORIES

Tel: 01759 301142 Fax: 01759 301143 Email: sales@wolflabs.co.uk
Website: www.wolflabs.co.uk

Test Tube Racks – Tube Peg Racks



Fibreglass-Reinforced Peg Racks Strong, tough, versatile, and attractive Autoclavable, fibreglass-reinforced polypropylene racks are ideal for bench work or storage of test tubes. Spacious peg design allows easy viewing of tube contents. User can also turn tubes upside-down and use as a drying rack! Fibreglass reinforcement keeps pegs stiff and straight through repeated autoclaving. Includes embossed indexing, ribbed grips, and drainage holes Pack of 2

Cat No	Tube Size	Wells	Colour
HS24311B	14-17mm	50	Blue
ZXZHS24312A	10-13mm	96	White
ZXZHS24312B	10-13mm	96	Blue