



Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax:01759 301143

sales@wolflabs.co.uk



If the literature below doesn't answer all your questions, click these links to:

[Request details on an alternative product of the same type](#)

[Request details on accessories](#)


[Request details on other product categories](#)



WOLF LABORATORIES

Tel: 01759 301142 Fax: 01759 301143 Email: sales@wolflabs.co.uk
Website: www.wolflabs.co.uk

Test Tube Racks – Wire 16mm



Hold up to 72 or 108 tubes

Chemical-resistant, epoxy-coated steel racks feature a slotted bottom and reinforced edges. Tubes sit securely in the square holes and the middle rung provides added support. Use smaller 72-well rack where space is at a premium. Autoclavable.

Cat No.	Tubes	Wells	Color	L x W x H, mm	Units
ZXZHS23108	16 mm	108	Red	241 x 188 x 65	Each
ZXZHS23072	16 mm	72	Blue	237 x 127 x 65	Each