

WOLF LABORATORIES

Tel: 01759 301142 Fax: 01759 301143 Email: sales@wolflabs.co.uk
Website: www.wolflabs.co.uk

Test Tube Wire Racks 15ml & 50ml



Sturdy racks for transport or storage

These chemical-resistant, epoxy-coated steel wire racks have a slotted bottom for better drainage and handling. Ideal for conical centrifuge tubes but can be used to hold other types and sizes of tubes. The 15 mL rack has 20 mm holes; 50 mL rack has 30 mm square holes
Autoclavable

Cat No.	Tubes	Wells	Colour	L x W x D mm
ZXZHS23215	15mL	15	Blue	127 x 76 x 76
ZXZHS232461	50mL	8	Blue	147 x 80 x 89