

# WOLF LABORATORIES

Tel: 01759 301142 Fax: 01759 301143 Email: sales@wolflabs.co.uk  
Website: www.wolflabs.co.uk

## Big Digit Dual Timer/Clock



Dual display permits you to simultaneously display a clock and timer, or two independent timers (up to 24 hours). Use one timer to count up, or use two to independently count down. Audible and visual alarm tells you "time is up". The large buttons and simple programming procedure make this instrument easy to use. Display tilts to allow freestanding use on your benchtop, or use the bracket (included) to mount on the wall next to your incubator. Includes one AAA battery. Measures 82 x 140 x 13mm

Cat No.	Description	Units
ZXZHS24780	Big Digit Dual Timer/Clock	Each