


WOLF LABORATORIES

Tel: 01759 301142 Fax: 01759 301143 Email: sales@wolflabs.co.uk
Website: www.wolflabs.co.uk

Tube Rack 20 Well



Interlock multiple racks to hold more tubes

Polypropylene racks lock together on each side. Easily remove tubes using finger slot on front of rack. Sturdy rack has internal ribs for stability. Racks have large molded coordinates for easy tube identification. Autoclavable. Measures 213 x 90 x 50 mm.
Pack of 6.

Part Number	Description	Units
ZXZHS4284C	20-Well Tube Racks, White	Pck of 6