

# WOLF LABORATORIES

Tel: 01759 301142 Fax: 01759 301143 Email: [sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)  
Website: [www.wolflabs.co.uk](http://www.wolflabs.co.uk)

## Tube Rack 3 Way



### Snap together for more tube storage

Polypropylene racks connect to one another on either end. Convenient finger slots are molded into the rack making them easy to pick up and move. Lowest level has 24 wells to hold individual 0.2 mL tubes or strips of 8 or 12 tubes. Middle tier has 14 wells to hold 0.5 mL tubes. Top tier has 12 wells to hold 1.5 mL tubes. Autoclavable. Available in natural or assorted colors: blue, green, purple, yellow, and orange. Rack measures 164 x 123 x 62 mm. Pack of 6

Cat No.	Colour	Units
ZXZHS4289A	Natural	Pck of 6