

WOLF LABORATORIES

Tel: 01759 301142 Fax: 01759 301143 Email: sales@wolflabs.co.uk
Website: www.wolflabs.co.uk

Weighing Boats



Standard Weighing Boats

Our weighing boats are a good addition to our white flexible boats. Both flexible boats have smooth surfaces to provide accurate pour-outs with minimal sample loss. These polystyrene boats have a shallow shape with a wide, flat bottom to resist tipping. Rounded corners simplify transfer. Biologically inert and resistant to dilute and weak acids, aqueous solutions, alcohols, and bases. May be used with microwave moisture analyzers. Resist temperatures to 200°F (93°C). Pack of 500.

These polystyrene boats are shallow with a wide, flat bottom to resist tipping. Rounded corners simplify transfer. May be used with microwave moisture analyzers

Cat No.	Colour	Size	L x W x D mm	Units
HS1420A	White	Small	46 x 46 x 8	Pck. of 500
ZXZHS1420B	White	Medium	85 x 85 x 24	Pck. of 500
ZXZHS1420C	White	Large	140 x 140 x 22	Pck. of 500
ZXZHS1421A	Blue	Small	46 x 46 x 8	Pck. of 500
ZXZHS1421B	Blue	Medium	85 x 85 x 24	Pck. of 500
ZXZHS1421C	Blue	Large	140 x 140 x 22	Pck. of 500