

BIG SQUID IKAMAG® White



Magnetic stirrer in five motifs.

- Diameter: 160 mm
- Stirring quantity (H₂O) up to 1.5 l
- Recyclable materials
- Very good chemical resistance due to glass top and synthetic bottom made of Hytrel®
- Electronically controlled motor with infinitely variable speed from 0 – 1,500 rpm

Technical Data			
Number of stirring positions [-]	1	Stirring quantity max. per stirring position (H ₂ O) [l]	1.5
Stirring quantity max. (H ₂ O) [l]	1.5	Motor rating input [W]	2
Motor rating output [W]	1	Speed range [1/min]	0 – 1500
Stirring bar length max. [mm]	30	Set-up plate material [-]	glass
Set-up plate dimensions [mm]	Ø 160	Dimensions (W x H x D) [mm]	180 x 48 x 180
Weight [kg]	1	Permissible ambient temperature [°C]	5 – 40
Permissible relative moisture [%]	80	Protection class according to DIN EN 60529 [-]	IP 54
Voltage [V]	100 – 240		