

RW 11 basic "Lab egg"



RW 11 basic "Lab egg" - 2830000

Price: **147** Ex Works - United Kingdom

The small-sized stirrer from IKA is available in 5 different designs with attractive colors. Since splashes come from below, glass, a material which is particularly resistant to chemicals, was used on this end of the stirrer.

The "Lab egg" is powerful: The power output of the agitator is designed for a maximum of 2 L of H₂O. Vane stirrer R 1001 and the extension arm are included as part of the package.



R 1001 spare paddle stirrer - 0527400

Price: **34.5** Ex Works - United Kingdom

Spare paddle stirrer for stirring motor RW 11 basic.

Stirrer diameter: 34 mm

Shaft diameter: 4 mm

Shaft length: 160 mm