

## Attachments for IKA KS260 Range

### AS 260.2 Fixing clip attachment

For shaking flasks, Erlenmeyer flasks and bottles with a round crosssection (without fixing clips).



Ident. No. **3115500**

Price **GBP 153,00**

### AS 260.3 Dish attachment

For smooth movement for cell cultures, nutrient media in Petri dishes, culture bottles and vessels with a low center of gravity. With integrated slip-resistant foil (PP).



Ident. No. **3120600**

Price **GBP 30,10**

### AS 260.5 Separating funnel attachment

For shaking out, salting out, extracting, eluting, enriching. The 3 clamping rolls (included in delivery) are height-adjustable for adaption to different separating funnel sizes. The separating funnels are secured with O-rings (5 O-rings included).



Ident. No. **3120900**

Price **GBP 171,00**

### AS 260.1 Universal attachment

For various types of vessels. Infinitely variable clamping rolls allow universal adaptation to the vessels.



Ident. No. **8017400**

Price **GBP 217,00**