

HS 260 basic



- Compact, flat shaker with ideal swivel motion, for a maximum shaking weight of 7.5 kg.
- Electronic adjustment of speed and timer
 - LED display for speed and time adjustment
 - Wide range of attachment combinations makes it possible to use almost all shapes and sizes of vessels
 - Attachments are not included in delivery, please order separately

Technical Data			
Type of movement [-]	reciprocating	Orbital diameter [mm]	20
Permissible shaking weight (incl. attachment) [kg]	7.5	Motor rating input [W]	45
Motor rating output [W]	10	Permissible ON time [%]	100
Speed range [1/min]	20 – 300	Speed display [-]	LED line
Timer [-]	yes	Timer display [-]	LED line
Time setting range [min]	5 – 50	Operating mode [-]	timer and continuous operation
Shaking table lock [-]	no	Dimensions (W x H x D) [mm]	360 x 100 x 420
Weight [kg]	8.1	Permissible ambient temperature [°C]	5 – 50
Permissible relative moisture [%]	100	Protection class according to DIN EN 60529 [-]	IP 21
RS 232 interface [-]	no	Analog output [-]	no
Voltage [V]	220 – 240	Frequency [Hz]	50/60
Power input [W]	45		