

RO 15 power IKAMAG®



Multi-position magnetic stirrer with 15 stirring positions, without heating. The stainless steel surface covers the unit allowing easy cleaning and providing protection against the penetration of liquids.

- Optimum use of laboratory space
- Including removable PUR cover

Technical Data			
Number of stirring positions [-]	15	Stirring quantity max. per stirring position (H ₂ O) [l]	0.4
Stirring position distance [mm]	90	Speed deviation of individual stirring positions [%]	5
Stirring quantity max. (H ₂ O) [l]	6	Motor rating input [W]	21.6
Motor rating output [W]	5.4	Speed range [1/min]	0 – 1100
Set-up plate material [-]	stainless steel 1.4301	Set-up plate dimensions [mm]	270 x 450
Dimensions (W x H x D) [mm]	272 x 65 x 552	Weight [kg]	4.7
Permissible ambient temperature [°C]	5 – 40	Permissible relative moisture [%]	80
Voltage [V]	230 / 115	Frequency [Hz]	50/60
Power input [W]	30		