

IKA Vortex Mixer



TTS 2 - 9008200

- Universal shaker for test tubes and small vessels
- Continuous operation or automatic mode of operation by pressure on the top attachment
- A large selection of attachments allows for an extensive range of applications
- Price incl. attachments 1.31 - 1.34

Shaking motion:	circular
Shaking orbit:	Ø 4,5 mm
Perm. agitation weight (with attachment):	0,5 kg
Motor rating input / output:	18 / 10 W
Permissible ON-time:	100 %
Speed range (infinitely adjustable):	0 - 2.500 rpm
Speed display:	Scale
Dimensions (W x D x H):	115 x 225 x 60 mm
Weight:	3,5 kg

Technical Data

TTS 2	
Max. Load With Attachment	.5 kg
Max. Shaking Frequency	2500 1/min
Shaking Motion	orbital rocking




MS 1.31 Test tube insert - L001840

Test tube insert for the test tube support MS 1.3. It can be used for 14 tubes with a diameter of 8 mm.



MS 1.32 Test tube insert - L001850

Test tube insert for the test tube support MS 1.3. It can be used for 6 tubes with a diameter of 10 mm.



MS 1.33 Test tube insert - L001860

Test tube insert for the test tube support MS 1.3. It can be used for 4 tubes with a diameter of 15 mm.



MS 1.34 Test tube insert - L001830

Test tube insert for the test tube support MS 1.3. The foam insert comes without bore holes.