

IKA DI 18 Homogeniser



DI 18 basic - L005118

- Attractively priced, narrow disperser
- A large speed range makes it possible to work at high circumferential speed
- Accepts plug-in dispersing tools

Speed range:	6.000 - 24.000 rpm
Motor capacity input/output:	300 / 160 W
Adjustment of speed:	infinitely variable
Speed display:	Scale
Dimensions (W x D x H):	78 x 78 x 220 mm
Weight:	1,2 kg

Dispersing tool, Plate stand and Boss head are not included und must be ordered separately, see Accessories.

Technical Data

DI 18 basic	
Max. Speed	24000 1/min
Max. Working Volume	0.1 L
Min. Batch Volume	.001 L
Process Type	batch

Accessories

S 18 N - 10 G - L004639



Dispersing tool for the T 18 Basic.
Working range: 1-100 ml



R 182 Boss head - 2657700

This boss head is used for the safe fixing of drives to the stands. It resists strongest mechanical demands. The new design has been awarded.

Stand diameter: 6-16 mm

Bracket diameter: 6-16 mm

***Recommended Accessory* for the R 1822 Plate stand.**



R 1826 Plate stand - 3160100