


# IKA® Mechanical accessories

## Fixing elements



**H 44 Boss head clamp**

Clamping range - stand	10 - 11 mm
Clamping range - extension arm	11 mm
Material	cast aluminium

Ident. No.  
**2437700**



**R 182 Boss head clamp**

Clamping range - stand	6 - 16 mm
Clamping range - extension arm	6 - 16 mm
Material	cast aluminium

Ident. No.  
**2657700**



**R 270 Boss head clamp**

Clamping range - stand	25 - 36 mm
Clamping range - extension arm	5 - 21 mm
Material	cast aluminium


Ident. No.  
**2657800**



**R 271 Boss head clamp**  
Specialized clamp with openings for the stands R 2722 and R 2723 as well as extensions with Ø 16 mm.

Clamping range - stand	34 mm
Clamping range - extension arm	16 mm
Material	cast aluminium


Ident. No.  
**2664000**



**RH 3 Strap clamp**  
For securing vessels against walls or for synchronized rotation during stirring or dispersing.

For stand diameter	8 - 16 mm
For vessel diameter	40 - 300 mm

Ident. No.  
**3008600**



**RH 5 Strap clamp**  
For securing vessels against walls or for synchronized rotation during stirring or dispersing, incl. boss head clamp R 270.

For stand diameter	25 - 36 mm
For vessel diameter	40 - 300 mm

Ident. No.  
**3159000**