

RH basic 2 IKAMAG®



RH basic 2
Economic magnetic stirrer with stainless steel heating plate.
- Fixed safety circuit 400 °C
- Soft-start stirring motor

Technical Data			
Number of stirring positions [-]	1	Stirring quantity max. per stirring position (H2O) [l]	10
Stirring quantity max. (H2O) [l]	10	Motor rating input [W]	15
Motor rating output [W]	2	Speed range [1/min]	100 - 2000
Heat output [W]	400	Heating temperature range [°C]	room temp. - 320
Heat control [-]	scale 1 - 6	Speed control [-]	scale 0 - 6
Fixed safety circuit [°C]	400	Set-up plate material [-]	stainless steel 1.4301
Set-up plate dimensions [mm]	Ø 125	Dimensions (W x H x D) [mm]	168 x 105 x 220
Weight [kg]	2.4	Permissible ambient temperature [°C]	5 - 40
Permissible relative moisture [%]	80	Protection class according to DIN EN 60529 [-]	IP 21
Voltage [V]	220 - 240 / 115 / 100	Frequency [Hz]	50/60
Power input [W]	415		