

## Data Sheet

Mixing / Magnetic Stirrers



## Topolino

Extremely convenient magnetic mini-stirrer for mixing quantities up to 250 ml.

- Durable, brushless motor
- Infinitely variable speed range
- High magnetic adhesion

Accessories: IKAFLON® 10 Magnetic stirring bar, IKAFLON® 15 Magnetic stirring bar, IKAFLON® 20 Magnetic stirring bar, IKAFLON® 25 Magnetic stirring bar, IKAFLON® 30 Magnetic stirring bar, TRIKA® 25 Magnetic stirring bar, RSE Stirring bar remover

### Technical Data

Number of stirring positions	1
Stirring quantity max. per stirring position (H <sub>2</sub> O) [l]	0,25
Stirring quantity max. (H <sub>2</sub> O) [l]	0,25
Motor rating input [W]	1
Motor rating output [W]	0,8
Speed range [1/min]	300 - 1800
Stirring bar length max. [mm]	30
Speed control	stepless
Set-up plate material	synthetic (PP)
Set-up plate dimensions [mm]	Ø 80
Dimensions (W x H x D) [mm]	95 x 37 x 115
Weight [kg]	0,32
Permissible ambient temperature [°C]	5 - 40
Permissible relative moisture [%]	80
Protection class according to DIN EN 60529	IP 21
RS 232 interface	no
Analog output	no
Voltage [V]	100 - 240
Frequency [Hz]	50/60
Power input [W]	10
<b>Ident. No.</b>	<b>3368000</b>