


## Multi-position magnetic stirrers with heating



### RT 5 power IKAMAG®

The RT 5 power is a high-performance multi-position magnetic stirrer with 5 stirring positions and integrated temperature control plate. Precise temperature distribution on the heating plate makes it possible to perform series experiments.


- Simultaneously operating stirrers
- Absolute consistency over sample conditions for the individual samples

#### Accessories (Page):

IKAFLO<sup>®</sup> Stirring bar (19),  
 TRIKA<sup>®</sup> Stirring bar (19),  
 RSE Stirring bar remover (19)

<b>Magnetic stirrer</b>	
Stirring positions	5
Max. stirring quantity per stirrer (H <sub>2</sub> O)	0.4 l
Distance between stirring places	90 mm
Motor rating input / output	7.2 W / 1,8 W
Speed display scale (1 - 10)	scale (1 - 10)
Speed range	0 - 1,100 rpm
Deviation for individual stirring positions	5 %
Max. magnetic bar (L x Ø)	30 x 8 mm
<b>Heating function</b>	
Heat output	175 W
Temperature range	RT - 120 °C
Heat control scale (1 - 10)	scale (1 - 10)
Temperature consistency in the medium	± 2 K
<b>Heating plate</b>	
Material	silicone
Dimensions	135 x 450 mm
<b>General data</b>	
Dimensions (W x D x H)	138 x 552 x 65 mm
Weight	3 kg
Permissible ambient temperature	5 - 40 °C
Permissible humidity	80 %
Protection class acc. to DIN EN 60529	IP 42

Ident. No.	Ident. No.
<b>2930300</b> 230 V 50/60 Hz	<b>2930301</b> 115 V 50/60 Hz



### RT 10 power IKAMAG®

Same features as RT 5 power, but with 10 stirring positions.

#### Accessories (Page):

IKAFLO<sup>®</sup> Stirring bars (19),  
 TRIKA<sup>®</sup> Stirring bars (19),  
 RSE Stirring bar remover (19)

<b>Magnetic stirrer</b>	
Stirring positions	10
Motor rating input / output	14.4 W / 3,6 W
<b>Heating function</b>	
Heat output	375 W
<b>Heating plate</b>	
Dimensions	195 x 450 mm
<b>General data</b>	
Dimensions (W x D x H)	198 x 552 x 65 mm
Weight	4.2 kg

Ident. No.	Ident. No.
<b>2930500</b> 230 V 50/60 Hz	<b>2930501</b> 115 V 50/60 Hz

### RT 15 power IKAMAG®

Same features as RT 5 power, but with 15 stirring positions


#### Accessories (Page):

IKAFLO<sup>®</sup> Stirring bars (19),  
 TRIKA<sup>®</sup> Stirring bars (19),  
 RSE Stirring bar remover (19)

<b>Magnetic stirrer</b>	
Stirring positions	15
Motor rating input / output	21.6 W / 5,4 W
<b>Heating function</b>	
Heat output	580 W
<b>Heating plate</b>	
Dimensions	285 x 450 mm
<b>General data</b>	
Dimensions (W x D x H)	288 x 552 x 65 mm
Weight	6 kg

Ident. No.	Ident. No.
<b>2930700</b> 230 V 50/60 Hz	<b>2930701</b> 115 V 50/60 Hz

## Multi-position magnetic stirrers without heating



### RO 5 power IKAMAG®

Multi-position magnetic stirrer with 5 stirring positions, without heating. The stainless steel surface covers the unit allowing easy cleaning and providing protection against the penetration of liquids.


- Optimum use of laboratory space
- Including removable PUR cover

#### Accessories (Page):

IKAFLO® Stirring bars (19),  
 TRIKA® Stirring bars (19),  
 RSE Stirring bar remover (19)

Magnetic stirrer	
Stirring positions	5
Max. stirring quantity per stirrer (H <sub>2</sub> O)	0.4 l
Distance between stirring places	90 mm
Motor rating	
input / output	7.2 W / 1.8 W
Speed display	scale (1 - 10)
Speed range	0 - 1,100 rpm
Deviation for individual stirring positions	5 %
Max. magnetic bar (L x Ø)	30 x 8 mm
Set-up plate	
Material	stainl. steel (1.4301)
Dimensions	120 x 450 mm
General data	
Dimensions (W x D x H)	122 x 552 x 65 mm
Weight	2.3 kg
Permissible ambient temperature	5 - 40 °C
Permissible humidity	80 %
Protection class acc. to DIN EN 60529	IP 42

Ident. No.	Ident. No.
<b>2930200</b> 230 V 50/60 Hz	<b>2930201</b> 115 V 50/60 Hz



### RO 10 power IKAMAG®


Same features as RO 5 power, but with 10 stirring positions.

#### Accessories (Page):

IKAFLO® Stirring bars (19),  
 TRIKA® Stirring bars (19),  
 RSE Stirring bar remover (19)

Magnetic stirrer	
Stirring positions	10
Motor rating	
input / output	14.4 W / 3.6 W
Set-up plate	
Dimensions	180 x 450 mm
General data	
Dimensions (W x D x H)	182 x 552 x 65 mm
Weight	3.2 kg

Ident. No.	Ident. No.
<b>2930400</b> 230 V 50/60 Hz	<b>2930401</b> 115 V 50/60 Hz



### RO 15 power IKAMAG®

Same features as RO 5 power, but with 15 stirring positions.

#### Accessories (Page):

IKAFLO® Stirring bars (19),  
 TRIKA® Stirring bars (19),  
 RSE Stirring bar remover (19)

Magnetic stirrer	
Stirring positions	15
Motor rating	
input / output	21.6 W / 5.4 W
Set-up plate	
Dimensions	270 x 450 mm
General data	
Dimensions (W x D x H)	272 x 552 x 65 mm
Weight	4.7 kg

Ident. No.	Ident. No.
<b>2930600</b> 230 V 50/60 Hz	<b>2930601</b> 115 V 50/60 Hz