



### H 30 Oil bath attachment

The oil bath attachments H 29 and H 30 can be used as oil baths together with an IKAMAG® magnetic stirrer with heating or with an IKATHERM® heating plate having a diameter of 135 mm.

- Positioning border prevents sliding on the heating plate
- Safety grip protects from bath attachment burns
- The bath attachment is made of aluminum. This ensures good heat transfer and quick heating-up of the tempering medium
- Easy cleaning
- The bath attachments can only be used as an oil bath
- Inner diameter 136-190 mm
- Height 110 mm
- Volume 1,5 Liter