


ilShin Up-Right Freezers and their exact dimensions




| Dimensions (mm) Up-Right models | | | | | |
|---------------------------------|--------|--------|--------|--------|--------|
| MODEL | DF8510 | DF8514 | DF8517 | DF8520 | DF8524 |
| Height 1 (top to bottom) | 1750 | 1979 | 1979 | 1979 | 1979 |
| Height 2 (without castors) | 1670 | 1955 | 1955 | 1955 | 1955 |
| Height 3 (Internal) | 960 | 1280 | 1280 | 1280 | 1280 |
| Height 4 (Between shelves) | 290 | 290 | 290 | 290 | 290 |
| Deep 1 (Front to Back) | 865 | 865 | 1015 | 1015 | 1120 |
| Deep 2 (without control panel) | 815 | 815 | 970 | 970 | 1080 |
| Deep 3 (without door) | 725 | 725 | 880 | 880 | 990 |
| Deep 4 (Internal) | 480 | 480 | 630 | 630 | 735 |
| Deep 5 (temperature probe) | 25 | 25 | 25 | 25 | 25 |
| Wide 1 (internal)* | 595* | 595* | 595* | 740* | 740* |
| Wide 2 (without door lock) | 890 | 890 | 890 | 1040 | 1040 |
| Wide 3 (external)** | 920 | 920 | 920 | 1070 | 1070 |

*The models with a stainless steel interior (DF8510S, DF8514S, DF8517S, DF8520S and DF8524S) have a reduced width of 12 mm.
(for example internal width of DF8510S is $595 - 12 = 583$ mm)

** DF8517 and DF8520 pass through a < 90 cm wide doorway with door open, front panels clicked off, and rolling the freezer from side to side.

ilShin Chest Freezers and their exact dimensions



| Dimensions (mm) Chest models | | | | | |
|------------------------------------|--------|--------|--------|--------|--------|
| MODEL | DF9007 | DF9010 | DF9014 | DF9017 | DF9020 |
| Height 1 (top to bottom) | 1010 | 1010 | 1010 | 1010 | 1010 |
| Height 2 (internal) | 665 | 665 | 665 | 665 | 665 |
| Height 3 (internal + extra frame)* | 695 | 695 | 695 | 695 | 695 |
| Deep 1 (Front to Back) | 790 | 790 | 790 | 790 | 790 |
| Deep 2 (without lid / hinges) | 740 | 740 | 740 | 740 | 740 |
| Deep 3 (internal) | 490 | 490 | 490 | 490 | 490 |
| Deep 4 (temperature probe) | 25 | 25 | 25 | 25 | 25 |
| Wide 1 (internal) | 560 | 860 | 1160 | 1460 | 1760 |
| Wide 2 (without door lock) | 1250 | 1550 | 1850 | 2150 | 2450 |
| Wide 3 (external) | 1280 | 1580 | 1880 | 2180 | 2480 |

*With use of a special ordered frame the internal height can be increased with 30 mm (for high storage racks)

| Model no: | Description: | Price: |
|-----------|---|----------|
| SF9007 | Stainles steel frame to increase chest height of DF9007 | € 215,00 |
| SF9010 | Stainles steel frame to increase chest height of DF9010 | € 215,00 |
| SF9014 | Stainles steel frame to increase chest height of DF9014 | € 215,00 |
| SF9017 | Stainles steel frame to increase chest height of DF9017 | € 215,00 |
| SF9020 | Stainles steel frame to increase chest height of DF9020 | € 215,00 |

Shipping Dimensions



| Model no: | Description: | Weight (kG): | W x D x H (mm): |
|--|-------------------------------------|--------------|--------------------|
| <i>Chest - 80 °C Freezers</i> | | | |
| DF9007 | - 80 CHEST FREEZER, CAPACITY-191L | 231 | 1680 x 1050 x 1160 |
| DF9010 | - 80 CHEST FREEZER, CAPACITY-295L | 256 | 1880 x 1050 x 1160 |
| DF9014 | - 80 CHEST FREEZER, CAPACITY-392L | 281 | 2100 x 1050 x 1160 |
| DF9017 | - 80 CHEST FREEZER, CAPACITY-492L | 306 | 2300 x 1050 x 1160 |
| DF9020 | - 80 CHEST FREEZER, CAPACITY-593L | 331 | 2740 x 1050 x 1160 |
| <i>Upright - 80 °C Freezers</i> | | | |
| DF8510 | -80 UP-RIGHT FREEZER, CAPACITY-276L | 314 | 1150 x 1050 x 1900 |
| DF8514 | -80 UP-RIGHT FREEZER, CAPACITY-369L | 354 | 1150 x 1050 x 2200 |
| DF8517 | -80 UP-RIGHT FREEZER, CAPACITY-484L | 400 | 1050 x 1150 x 2200 |
| DF8520 | -80 UP-RIGHT FREEZER, CAPACITY-605L | 430 | 1155 x 1150 x 2200 |
| DF8524 | -80 UP-RIGHT FREEZER, CAPACITY-730L | 470 | 1200 x 1150 x 2200 |
| <i>Flake Ice Maker</i> | | | |
| IF300-150 | FLAKE ICEMAKER | 166 | 970 x 970 x 1210 |

ilShin[®] Europe

Kryptonstraat 11
 6718 WR EDE
 The Netherlands
 Tel. +31 318 64 64 81
 Fax.+31 318 646034
www.ilshineurope.com