

VXS 80 - Personal -86°C Freezer

› VXS 80
Personal ULT

Our Jouan VXS 80 personal-sized freezer fits under or on top of your benchtop, ensuring quick direct access to your samples. Thermo offers a choice of configurations to optimize storage. Removable shelves and racks of different heights allow you to make the best use of space.

Product Features

- Space-efficient size keeps your samples within easy reach
- Freezer is quiet-running and dissipates little heat to ensure comfort in the lab
- High-performance microprocessor control, high temperature alarm, diagnostic alarm menu and key lock are standard
- Optimum cooling power provides rapid freezing and quick temperature recovery after door closings
- Optional CO₂ back-up is available



Customer Services
Phone 0870 609 9203
Fax 0870 609 9202
sales.labequip.uk@thermo.com
www.thermo.com



Easy Access and Optimized Use of Space



Range of Racks

Application

The Jouan VXS 80 personal freezer is

ideal for small cleanrooms and all culture laboratories; samples can be

close to the bench, easily within reach without leaving the room.

Specifications

Model No	VXS 80
Capacity cu ft (ltrs)	2.8 (80)
No of Cryoboxes (50 / 70 mm)	54 / 45
Temperature Range	-65°C to -86°C
External Dimensions H x W x D	33.5 x 25.6 x 26.8 (in) 850 x 650 x 680 (mm)
Internal Dimensions H x W x D	15.6 x 17.9 x 17.9 (in) 395 x 455 x 455 (mm)
Weight lbs (kg)	243 (110)

Ordering Information

Cat No	Description	Power Requirements
51501200	VXS 80 Personal -86°C Freezer	230V, 50 Hz

Accessories

Cat No	Description
51500599	CO ₂ safety back-up system
51500795	7 day chart recorder (independent - battery operated)
88340801	Extra shelf kit

Software

Cat No	Description
51501300	Freezer Monitoring System (FMS) software - 230V
51501301	Freezer Monitoring System (FMS) software - 208-220V

Consumables

Cat No	Description
51500518	Set of 5 pens for 7-day recorder
51500519	Set of 100 discs for 7-day recorder