

## Walk-In Fume Cupboard

Designed to accommodate tall test assemblies and to provide high headroom in full hazard and protection environments.

Features include a double sash assembly designed to provide independent sash operation and services to match the user's requirements. Options include a hinged or removable intermediate work surface.

Combining functionality with design, the walk-in is available in a range of colours to co-ordinate with surrounding furniture.



## Bench Mounted Fume Cupboard

Providing full hazard and containment protection using proven designs, these units can be bench mounted on existing furniture or supplied with an integral steel support unit and cupboards.

The cupboards are available in a range of colours to match laboratory furniture and worktops, with services and linings specified to meet your requirements.



## **Step-In Fume Cupboard**

Benefitting from a similar design to the walk-in fume cupboard but incorporating a low level work surface 400mm above floor level to allow for sinks and drip cups.

Step-in Fume Cupboards are frequently specified for use with distillation equipment. Services and linings are available to user's specifications. The unit is available in a range of colours to complement surrounding furniture.



## **Graduate Fume Cupboard**

A cost effective unit for use in low hazard laboratories such as schools and colleges, which complies with the current D of E 'Fume Cupboards for Schools - Bulletin 88.'

The Graduate is available for bench mounting onto existing furniture, or free standing - fitted with an integral steel support unit.

