

**BLOODBANK.** The Lec Medical bloodbank range complies with British standard BS 4376-1:1991 and offers a variety of options from the more compact sizes (130 to 250 litres) to much larger models (700 to 1500 litres).

**Compact Sized Features**

Microprocessor with digital display giving accurate temperature readings.

High/low temperature alarm system.

Visual indication of alarm occurrence after the event.

48 hour battery back-up in the event of a power failure.

Internal ventilation.

Internal light.

7 day chart recorder with 1.5 volt battery circuit to guard against power failure.

Supplied with one year's supply of charts and two pen nibs.

**Larger Sized Features**

Electric heaters evaporate condensed water to prevent flooding.

Microprocessor with digital display giving accurate temperature readings.

High/low temperature alarm system.

Electric heaters for automatic defrost.

Visual indication of alarm occurrence after the event.

48 hour battery back-up in the event of a power failure.

Internal ventilation.

Internal light.

Adjustable feet.

7 day chart recorder with 1.5 volt battery circuit to guard against power failure.

Supplied with one year's supply of charts and two pen nibs.

All bloodbanks are made to order and delivery will usually take between 6 to 8 weeks.

**BB36**

SKU 444447541  
 Product dimensions (mm) 600w x 700d x 1250h  
 Product weight (kg) 92  
 Gross Capacity (litres) 130 (32-40 packs)  
 Drawers 2

**BB72**

SKU 444447605  
 Product dimensions (mm) 600w x 700d x 1810h  
 Product weight (kg) 122  
 Gross Capacity (litres) 250 (65-79 packs)  
 Drawers 4

**BB90**

SKU 444448250  
 Product dimensions (mm) 600w x 700d x 1810h  
 Product weight (kg) 130  
 Gross Capacity (litres) 250 (81-99 packs)  
 Drawers 5

**BB264**

SKU 444448249  
 Product dimensions (mm) 740w x 815d x 1998h  
 Product weight (kg) 185  
 Gross Capacity (litres) 700 (240-290 packs)  
 Drawers 6

**BB352**

SKU 444447603  
 Product dimensions (mm) 740w x 815d x 1998h  
 Product weight (kg) 195  
 Gross Capacity (litres) 700 (320-390 packs)  
 Drawers 8

**BB528**

SKU 444447604  
 Product dimensions (mm) 1480w x 815d x 2020h  
 Product weight (kg) 290  
 Gross Capacity (litres) 1500 (530-580 packs)  
 Drawers 12

BB36



BB72



BB90



BB264



BB352



BB528

