

**WARD.** Lec Medical's range of ward refrigeration offers something for all requirements. From the compact WR207 for smaller wards through to the larger WR907 which is suitable for even the busiest of wards.

#### Control and Efficiency

The neat, integral digital control panel on the outside of the fridge is easy to use and makes checking the internal temperature straightforward. As there is no need to open the door to check the temperature, it keeps the internal air temperature more stable.

All three ward fridges have fans delivering optimum temperature stability and a faster temperature recovery after door openings.

Automatic defrost ensures more efficient cooling.

#### Security

Security is essential on all wards so all Lec Medical fridges come complete with door locks and two keys ensuring that only authorised individuals are able to access them.

A temperature alarm, which sounds if the temperature deviates outside of the set temperatures, comes as standard, reducing the chance of ruining the contents of the refrigerator.

The auxiliary output\*, gives the user the option of connecting support systems such as a mains failure alarm signal.

#### Flexibility

Reversible doors across all models give greater flexibility where space is at a premium or where positioning the product is specific.

Adjustable feet are also important when positioning the product on uneven surfaces.

Removable shelves and door furniture\* gives maximum storage capacity.

#### Environmentally Friendly

Low energy consumption.

CFC and HFC free.

#### WR207

SKU	444448570
Product dimensions (mm)	502w x 540d x 660h
Packed dimensions (mm)	520w x 595d x 675h
Product weight (kg)	27.5
Packed weight (kg)	30.5
Net capacity (litres)	67
Gross capacity (litres)	82

#### WR507

SKU	444448524
Product dimensions (mm)	558w x 600d x 863h
Packed dimensions (mm)	600w x 610d x 900h
Product weight (kg)	34
Packed weight (kg)	35.5
Net capacity (litres)	137
Gross capacity (litres)	153

#### WR907

SKU	444448525
Product dimensions (mm)	595w x 600d x 1550h
Packed dimensions (mm)	625w x 700d x 1610h
Product weight (kg)	81.5
Packed weight (kg)	84
Net capacity (litres)	275
Gross capacity (litres)	286

WR207



WR507



WR907

