


## LER 500-40 Laboratory Freezer ( -20°C to -40°C)



| Dimensions    | Height | Width | Depth | Weight | Capacity  |
|---------------|--------|-------|-------|--------|-----------|
| External (mm) | 1702   | 770   | 720   | 94Kg   |           |
| Internal (mm) | 1469   | 650   | 555   |        | 525 litre |

- Operating Temperature -20°C to -40°C
- Door Gasket easily replaced
- Door lock fitted
- Full length door handle
- 6 Static freezer shelves
- 100% CFC Free
- Optional Castor Base
- Exterior & Interior finish white ( Stainless Steel exterior option)
- Digital display of temperature and alarm condition with audible alarm.
- External volt free alarm contacts.
- Control keys lockable to prevent unauthorised adjustments.

### Options

- Temperature Recorder (Factory fitted)
- Inner doors
- Blood Plasma version generally conforming to BS 4376 P11
- Castor Base
- Stainless Steel exterior
- Storage boxes