



RAFR05202



RAFR44261



RAFR21261

Microprocessor control, together with ducted air flow and advanced auto defrost system provide class leading temperature stability and uniformity.

Options available (dependant on model): Glass door(s); stainless steel construction; Sample storage racking; Remote alarm contacts; Access ports; Internal power sockets with safety trip; Over and under temperature safety trips; Temperature recorders/loggers; Stacking kit for 129 litre models.

Advanced models may also be customised for specific applications – please ask us.

FREEZERS –20°C

Part Number	RAFR05202	RAFR05202GD	RAFR14261	RAFR21261	RAFG21261	RAFR44261	RAFG44261
Configuration	Underbench, solid door	Underbench, glass door	Free-standing solid door	Free-standing solid door	Free-standing glass door	Free-standing double solid door	Free-standing double glass door
Capacity (litres)	129	129	385	577	577	1245	1245
Exterior size (mm)							
High	850	850	1995	1995	1995	1995	1995
Wide	630	630	600	730	730	1460	1460
Deep	567	567	700	800	800	800	800
Sparkfree interior	X	X	X	X	X	X	X
Internal fan	X	X	✓	✓	✓	✓	✓
Door lock	✓	✓	✓	✓	✓	✓	✓
Temp display	✓	✓	✓	✓	✓	✓	✓
Alarm	✓	✓	✓	✓	✓	✓	✓
Auto defrost	X	X	✓	✓	✓	✓	✓
No. shelves	3	3	4	4	4	8	8
Operating temp (nominal)	-20°C	-20°C	-20°C	-20°C	-20°C	-20°C	-20°C

✓ – standard ○ – available as an optional accessory X – not available on this model