

# ADVANCED FRIDGES AND FREEZERS



RAFG21041



RAFR05042



RAFG44041

- **Precise microprocessor temperature control**
- **Digital temperature display**
- **Audible and visual temperature alarm**
- **Optimised air circulation and temperature uniformity**
- **Automatic defrosting with minimal temperature deviation**
- **Lockable castors on free-standing models**
- **Heavy duty all-metal construction**
- **Lockable doors with choice of glass or solid construction**
- **Adjustable heavy duty shelving**

**The Labcold Advanced** range of refrigerators and freezers, along with their Pharmacy counterparts, represent the highest levels of performance, constructional quality, robustness, and appearance available in the market place. Heavy duty all-metal construction, with standard or optional stainless steel (dependant upon model) provides strength and durability with ease of cleaning or sanitisation.

## REFRIGERATORS +4°C

Part Number	RAFR05042	RAFR05042GD	RAFR14041	RAFR21041	RAFG21041	RAFR44041	RAFG44041
<b>Configuration</b>	Underbench, solid door	Underbench, glass door	Free-standing solid door	Free-standing solid door	Free-standing glass door	Free-standing double solid door	Free-standing double glass door
<b>Capacity (litres)</b>	129	129	385	577	577	1245	1245
<b>Exterior size (mm)</b>							
<b>High</b>	850	850	1995	1995	1995	1995	1995
<b>Wide</b>	630	630	600	730	730	1460	1460
<b>Deep</b>	567	567	700	800	800	800	800
<b>Sparkfree interior</b>	X	X	X	X	X	X	X
<b>Internal fan</b>	✓	✓	✓	✓	✓	✓	✓
<b>Door lock</b>	✓	✓	✓	✓	✓	✓	✓
<b>Temp display</b>	✓	✓	✓	✓	✓	✓	✓
<b>Alarm</b>	✓	✓	✓	✓	✓	✓	✓
<b>Auto defrost</b>	✓	✓	✓	✓	✓	✓	✓
<b>No. shelves</b>	3	3	4	4	4	8	8
<b>Operating temp (nominal)</b>	+4°C	+4°C	+4°C	+4°C	+4°C	+4°C	+4°C

✓ - standard    ○ - available as an optional accessory    X - not available on this model