

COMBINATION FRIDGE/FREEZERS



RLLR05042



RLFF13245



RAFF18241



RLFF10243

- **Four models to suit most applications**
- **Sparkfree or fan circulated models**
- **Refrigerator sections automatic defrosting**
- **Space saving combined cabinet design**
- **Low energy models**
- **Operating temperatures +4/-20°C**
- **CFC and HCFC free**

These models provide the same features and performance as the Sparkfree, Economy or Advanced ranges from which they derive, but with the convenience of a single cabinet. Both refrigerated and frozen samples can be stored and accessed within the same small floor area.

Options available (dependant on model): lockable models; temperature alarm/display; sample storage racking.

Part Number	RLLR05042	RLFF10243	RLFF13245	RLFF13245/LK	RAFF18241
Configuration	Underbench	Free-standing	Free-standing	Free-standing, lockable	Free-standing
Capacity (litres)*	106/17	160/126	170/160	170/160	250/250
Exterior size (mm)					
High	850	1650	1985	1985	1995
Wide	550	550	600	625	730
Deep	620	610	630	630	800
Sparkfree interior	✓	✗	✓	✓	✗
Internal fan	✗	✗	✗	✗	✓
Door lock	✓	✗	✗	✓	✓
Temp display	○	✗	○	○	✓
Alarm	○	✗	○	○	✓
Auto defrost*	✓/✗	✓/✗	✓/✗	✓/✗	✓/✓
Low energy	✓	✓	✓	✓	✗
No. shelves/baskets*	2/0	3/4	4/5	4/5	2/2
Operating temp (nominal)*	+4/-20°C	+4/-20°C	+4/-20°C	+4/-20°C	+4/-20°C

✓ – standard ○ – available as an optional accessory ✗ – not available on this model
 * refrigerator section/freezer section