

ECONOMY FRIDGES AND FREEZERS



RLFR1804



RLVL03202



RLVF1825

- **Low cost laboratory units**
- **Underbench or free-standing models**
- **Auto defrosting refrigerators**
- **Operator controlled defrosting on freezers**

These models provide basic refrigerated storage at a lower price while still providing the confidence of the Labcold name and support. Lower cost production methods have been used, and certain convenience features omitted, in order to offer the best price possible without significant sacrifices in the important areas of function and reliability.

Accessories are not available with these products.

Part Number	RLPL04042	RLFR1804	RLVL03202	RLVF1825
Configuration	Under-bench Refrigerator	Free-standing Refrigerator	Under-bench Freezer	Free-standing Freezer
Capacity (litres)	117	505	90	505
Exterior size (mm)				
High	855	1710	855	1710
Wide	505	765	505	765
Deep	500	725	500	725
Sparkfree interior	X	X	X	X
Internal fan	X	✓	X	X
Door lock	X	✓	X	✓
Temp display	X	✓	X	✓
Alarm	X	X	X	X
Auto defrost	✓	✓	X	X
No. shelves/drawers	4	4	4	6
Operating temp (nominal)	+4°C	+4°C	-20°C	-20°C

✓ – standard X – not available on this model