

PRODUCT INFORMATION

Model: RECW RANGE

Type: Labcold Electronic Chart Recorder

Labcold chart recorders provide permanent records of storage temperatures in a wide variety of cabinets. Recorders have a PT1000 temperature sensor and are supplied suitable for wall mounting (fitting not included).



Recorder shown is panel mount version

Temperature ranges:	RECW0050:	0°C to +50°C
	RECW1020:	-10°C to +20°C
	RECW5000:	-50°C to 0°C
	RECW9010:	-90°C to +10°C

Sensor Dimensions:	50 x 6mm diameter
Power:	240V
Battery back up:	Fitted as standard
Pen Life:	Approx. 6 months
Exterior dimensions (mm):	144 x 144 x 67
Accuracy:	±1%
Recording period:	7 days per chart

Order reference: RECW****